Lovely Light Bleaching Paste is a gentle bleaching crème enriched with Jojoba natural oil, which is moisturizing and protective. By using the bleaching paste you have no formation volatile, irritant powders.

Preparation: Before usage, wear protective gloves up to the end of the service. Protect the hairline with a barrier crème, ie: It's Color Neutro can be used in this manner. Always mix in a non-metallic plastic or glass bowl.

Mixing Ratio: 1:2

Ex: 30 grams Lovely Light bleaching paste + 60ml Lovely Light Cream developer 10, 20 or 30 Volume. Mix with a brush until soft and creamy with no lumps. To adjust the mixture to the various types of hair for bleaching refer to the mixing chart.

Application: Apply the mixture on dry and always unwashed hair starting from parts which are most likely more difficult to lighten. If desired result requires an application on the scalp, never use a Lovely Light developer which is higher than 20 vol. Be sure that the scalp doesn't have irritations or other visible damage and verify the absence of hypersensitivity by performing a preliminary test.

Processing Time: Leave on the time required to achieve the desired result by often controlling the level of bleaching achieved. The final result will depend also on other factors such as: natural or cosmetic color, nature and structure of hair, quantities and volumes of developer used in mixture. Do not use above 104° F.

Rinsing: At the end of processing time, rinse thoroughly and shampoo with Easy Care Magical Color Shampoo. Follow up with a Magical Color Conditioner treatment.

Correction of warm reflexes: If necessary follow with a tonalization service with YouUp2 using 10.01 (10NA) or Level 11 Lightest Silver with a 1:2 mixing ratio with specific YouUp2 Energizer.



- On scalp bleach in a paste formula
- NO DUSTING so it's easier to mix
- Multiple application
- tube.

## **Table of Mixing**

Structure and level of hair to bleach	Lovely Light Dev.	Mixing Ratio	Processing Time	Performance of Highlight
Natural thick, difficult, dark natural or col- ored up to Lvl 4	10 vol	1:1	15 min	1-2 levels
	20 vol	1:1	30 min	2-3 levels
	30 vol	1:1	45 min	4-5 levels
Natural normal or thin, medi- um/ light natu- ral or colored from Lvl 4-7	10 vol	1:2	15 min	2-3 levels
	20 vol	1:2	30 min	4-5 levels
	30 vol	1:2	45 min	4-5 levels
Regrowth bleaching from Lvl 4-7	10 vol	1:2	15 min	2-3 levels
	20 vol	1:2	30 min	4-5 levels
Regrowth bleaching from Lvl 4-7	10 vol	1:2	15 min	3-4 levels
	20 vol	1:2	30 min	4-6 levels

5.1 oz tube